Description

Evetric® Bond is a light-curing, single-component dental adhesive for enamel and dentin in a bottle delivery form for repeated dispensing.

The adhesive is used:

- in the "self-etch" technique;
- with "selective enamel etching".

Composition

Evetric Bond contains methacrylates, ethanol, water, highly dispersed silicon dioxide, initiators and stabilizers.

Indications

- Direct, light-curing composite and compomer restorations.
- Direct core build-ups with light-, self- and dual-curing composites.
- Repair of fractured composite and compomer restorations.
- Desensitization of sensitive tooth necks.

Contraindication

- Do not use Evetric Bond if a patient is known to be allergic to any of the materials' ingredients or if the stipulated working technique cannot be used.
- Applications for which sufficient illumination is not ensured.
- Direct pulp capping.

Side effects

In rare cases, components of Evetric Bond may lead to sensitization. In these cases, Evetric Bond should no longer be used.

Interactions

- Phenolic substances may inhibit the polymerization of Evetric Bond (e.g. eugenol, salicyl derivates). Such materials should not be used ≤ 12 h prior to the application of Evetric Bond.
- Iron- and aluminium-based styptic agents must not be used shortly before or during the application of Evetric Bond.
- Bleaching agents or oxidizing disinfectants must not be used ≤ 7 days prior to the application of Evetric Bond.
- Large-dose fluoridation agents must not be used ≤ 7 days prior to the application of Evetric Bond.

Procedure**1. Isolation**

Ensure that the operating field is dry. Ideally, place a rubber dam (e.g. OptraDam® Plus).

2. Pulp protection / Base

In very deep, close-to-the-pulp cavities, selectively coat this area with a calcium hydroxide liner (e.g. ApexCal®) and subsequently cover it with pressure-resistant cement (glass ionomer cement, e.g. Vivaglass® Liner).

3. Conditioning with phosphoric acid gel for selective enamel etching (optional)

The bond to enamel can be additionally improved by selectively etching the enamel. Please observe the Instructions for Use of the phosphoric acid gel.

Apply phosphoric acid gel (e.g. Total Etch) to the enamel surface and allow to react for 15 to 30 seconds. Then thoroughly rinse off the gel with a strong jet of water for at least 5 seconds and dry with compressed air until the etched enamel surfaces appear chalky white.

4. Application of the adhesive

- Dispense the required quantity of Evetric Bond on a spot plate and gather it from there with a single-use applicator (e.g. Vivadent® Applicator Brush). Protect the adhesive from light (e.g. VivaPad®).
- Carefully close the bottle immediately after use.
- Entirely wet the tooth surfaces to be treated with Evetric Bond starting with the enamel.
- Rub the adhesive into the tooth surface to be treated for at least 20 seconds. Do not shorten this time. Merely spreading the adhesive on the tooth surface is not sufficient.
- Disperse Evetric Bond with oil- and moisture-free compressed air until a glossy, immovable film has developed.
- Light-cure Evetric Bond for 10 seconds using a light intensity of ≥ 500 mW/cm² (e.g. Bluephase® Style).
- Use a new disposable applicator for each new application.
- Refer to the Instructions for Use of the restorative materials for the further procedure.

5. Intraoral repair of fractured composite restorations with light-curing materials

- Roughen the surface to be repaired with a diamond finisher, clean with water spray and dry with oil- and moisture-free compressed air.
- Continue the further procedure as described in Point 4 (Application of the adhesive).

Notes

- Contamination at the bottle may be removed by e.g. wiping it with alcohol or another non-oxidizing disinfectant.
- Do not use the bottle in the mouth of the patient. Use a suitable application aid (e.g. VivaPad®) for reasons of hygiene.
- Unprepared enamel surfaces must be conditioned with phosphoric acid (e.g. Total Etch).

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Warning

- Evetric Bond is irritating. Avoid contact with the skin, mucous membrane and eyes.
- In case of skin contact, immediately wash with soap and copious amounts of water.
- If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician / ophthalmologist.
- In rare cases, contact with the skin may lead to sensitization towards the ingredients.
- Customary medical gloves do not provide protection against the sensitizing effect of methacrylates.

Notes regarding phosphoric acid

Phosphoric acid gel is etching. Avoid contact with the skin, mucous membrane and eyes (goggles are recommended for both patients and practitioners). If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician / ophthalmologist.

Storage

- Expiration date: see imprint on packaging.
- Immediately close the bottle with the original cap after use.
- Storage temperature: 2–28 °C / 36–82 °F.

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

The Material Safety Data Sheet is available online at www.ivoclarvivadent.com

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Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.